Growth Curve Report

$Randy\ L\ Coryell$ June 22, 2016

Contents

Overview of Growth Study
Strains:
Installation of R packages
Raw Data
Tidy the data
Overview of Growth Study
This report will highlight data aquired during a 8h growth study to "investigate the neutrality of Tn^2 antibiotic (CAM) neutral marker insertion in $Vibrio\ fischeri$
Strains:
ES114
ET401
EM17
Note All cultures were derived from a single colony, isolated on CAM plates (2.5ug/mL) and grown in liquid SWT (100mL) overnight. These overnight cultures were used to innoculate fresh SWT (100mL) media and grown at 225rpm and 28C for 1.5h. After taking a turbity measurement via spectrophotometry, each strain was diluted by placing a volume of this 1.5h culture into 100ml of fresh SWT. An initial OD was taken and the cultures were then placed into a shaking incubator @ 28C and 225RPM for the study.
I want to walk down the street and have people say 'There goes Roy Hobbs, the best to ever play the game'.
Robert Redford, The Natural

Installation of R packages

Certain R packages, available on the CRAN repository, are required for this analysis.

install.packages(c('reshape2', 'dplyr', 'ggplot2'))

Raw Data

In order to generate growth curves from the raw data, it first must be recorded as a CSV file in MSEXCEL and then loaded into memory in R.

```
rawdata <- read.csv("/Users/randycoryell/Growth Curve/Growth Curve/06212016_CAM_VS_NOCAMcurve.csv")
print(rawdata)</pre>
```

```
Time Temp ES114CAM ES114NO ET401CAM ET401NO EM17CAM EM17NO
##
## 1
        0
            28
                   0.016
                           0.062
                                     0.013
                                              0.060
                                                      0.082
                                                             0.090
## 2
       30
            28
                   0.257
                           0.258
                                     0.187
                                             0.257
                                                             0.126
                                                      0.146
## 3
       60
            28
                   0.389
                           0.430
                                     0.326
                                             0.355
                                                      0.143
                                                             0.135
                                     0.485
## 4
       90
            28
                   0.476
                           0.507
                                             0.513
                                                      0.143
                                                             0.134
## 5
      120
            28
                   0.618
                           0.671
                                     0.734
                                             0.747
                                                      0.151
                                                             0.130
## 6
      180
                   0.973
                                              1.230
            28
                           1.040
                                     1.120
                                                      0.106
                                                             0.131
## 7
      240
            28
                   1.450
                           1.530
                                     1.660
                                             1.560
                                                      0.175
                                                             0.210
## 8
      300
            28
                   1.730
                           1.880
                                     1.900
                                              1.860
                                                      0.164 0.436
## 9
      360
            28
                   2.170
                           2.310
                                     2.020
                                             2.000
                                                      0.178 0.795
```

Tidy the data

A reshaping of the data matrix is required in order to futher analyze the data set.

```
library(reshape2)
reshaped <- melt(rawdata, id=c("Time", "Temp"), variable.name="Strain",
value.name="OD600")</pre>
```

Summary statistics can be called upon for the reshaped data

summary(reshaped)

```
##
                                                    OD600
         Time
                          Temp
                                        Strain
              0.0
                                  ES114CAM:9
                                                        :0.0130
##
    Min.
           :
                    Min.
                            :28
                                                Min.
   1st Qu.: 60.0
##
                    1st Qu.:28
                                  ES114NO:9
                                                1st Qu.:0.1437
  Median :120.0
                    Median:28
                                  ET401CAM:9
                                                Median : 0.4095
## Mean
           :153.3
                    Mean
                            :28
                                  ET401NO :9
                                                Mean
                                                        :0.6902
##
    3rd Qu.:240.0
                    3rd Qu.:28
                                  EM17CAM :9
                                                3rd Qu.:1.1000
           :360.0
                                  EM17NO :9
## Max.
                    Max.
                            :28
                                                Max.
                                                        :2.3100
```

Now we need to take a quick look at the reshaped matrix to check for farmatting errors

head(reshaped)

```
Strain OD600
##
     Time Temp
## 1
            28 ES114CAM 0.016
        0
## 2
       30
            28 ES114CAM 0.257
## 3
       60
            28 ES114CAM 0.389
## 4
       90
            28 ES114CAM 0.476
      120
            28 ES114CAM 0.618
## 5
      180
            28 ES114CAM 0.973
## 6
```

When the data is replicated, there would need to be a nesting or grouping performed, this step would be useful for replicates in experimental evolution studies

```
grouped <- group_by(reshaped, Strain)</pre>
```

Sice these are single reads from the spectrophotometer, no confidence interval can be established, however an approximate growth curve can be produced.

library(ggplot2)

Warning: package 'ggplot2' was built under R version 3.2.4

print(ggplot(data=reshaped, aes(x=Time/60, y=OD600, color=Strain)) + geom_line() + labs(x="Time (Hours)

